



Notes


✨ Spark Testing ✨

Amazing things that worked wonders ❤️

- Programmer reads already created parameters via lead group amazingly
- Maïke is super happy about the import options
- Visually, the 3D viewer is really nice - chef's kiss

Suggested Changes 🐶

- Instead of erroring out, give a warning when clicking spark without patients in lead dbs
 - allow to erase added session withing the UI
 - would be good to give a warning when going back without saving (there's no feedback when saving or going back without saving clinical scores)
 - OSS-DBS settings button doesn't work when clicking
 - Hz wrote with lowercase
 - Maybe don't change the size of the UPDRS-III visualizations when setting clinical variables when hovering over them, maybe only increase the size when being clicked on?
 - When opening a different app (only via the desktop doc not switching between opened apps) and then going back to electron, javascript error appears
-



- Shouldn't the maximum amplitude defined in the left column match the amplitudes on the contacts on the electrode?

Unable to load the image

Upload Image

Errors

- Error in template excel for stimulation parameters (gives error and then is blank)



Alert

We found a problem with some content
in
'Stimulation_Parameters_Template.xlsx'.
Do you want us to try to recover as much
as we can? If you trust the source of this
workbook, click Yes.

Yes

No

Bugs 🤖

- Opening saved session → add stimulation parameters twice and all is gone (the first time it is clicked, something should pop, otherwise user will click the second time and spark crashes)
- Adding new clinical outcomes is not intuitive. When I click add it doesn't show whether it was added or not, but if you close it then it's showing as added

Create Clinical Score Type



ZscoreUPDRS

'Happy','Sad','Drunk'

separated by commas (Ex: 'Obsessions', 'Compulsions')

Close

Add